

PURIKER OS


Experimental shaker that can be used for a long time without failure due to its durability

**STANDARD SPECIFICATION**

MODEL	PURIKER OS
Speed range	30 ~ 300 rpm
Shaking Throw	20mm
Motor	BLDC Motor
Speed Control	Analog Jog Dial
Platform size (WxDxH)mm	300 x 300
Size (WxDxH)mm	320 x 420 x 120
Power Consumption	30W
Power	AC100 ~ 240V, 50/60Hz
Weight(Kg)	12.1
Cat. No.	CS-PKR-OS-Set

MAIN FEATURE

- Guaranteed long-term operation by applying BLDC motor
- Accurate speed control through feedback control system
- Provides excellent shaking quality through smooth operation
- Stable use even at high speed with low-height design
- Provide "auto stop" and "alarm function" to prevent excessive force from being applied to the product
- Provides "silicone pad" to prevent slipping

**ORDER INFORMATION**

DIVISION	ITEM	PURIKER OS
		Cat. No.
PRODUCT	Main Set	CS-PKR-OS-Set