


PURIMIX

Compact vortex mixer for simple mixing of liquid sample

**STANDARD SPECIFICATION**

MODEL		PURIMIX
Size	Out(W x D x H)mm	132 x 185 x 130
	Weight(Kg)	2.3
Mixing	Motion	Orbital
Material		Powder Coated Aluminum
Speed Range		0 to 3,300 RPM
Speed Control		Stepless
Operating Mode		Touch or Continuous
Throw		Circular Orbit, 4.00mm diameter
Driving Motor		Brushless Shade Motor
Shaking Load Capacity		0.5Kg
Permissible Temp. & Humidity		5°C~40°C, < 80%
Power		AC 230V±10%, 50/60Hz
Power Consumption		30W
Application Container/Number		Hand holding; any micro tube or conical tube
Cat. No.		CS-PMX-Set

FEATURE PLUS

3inch Rubber Platform Head



Pop-Off Cup Head




Remove Applicator

MAIN FEATURE



Multifunction Switch

Intuitive Power switch to control among three modes such as 'Continuous Action', 'Touch On' and 'Power Off'



Wide Application

Concavo-convex rubber head allows for mixing a variety of tubes with different shapes, sizes and materials.



Powerful performance and reliability

Powerful motor and robust steel case enable stable mixing.



Safety Assurance

Over heat and over current are not allowed to affect on a motor.



Changeable mixing speed

Mixing speed adjustment from slow speed to high speed (Max. 3,300rpm)



Long-life motor

Brushless shade motor has a long life cycle.

ORDER INFORMATION

DIVISION	ITEM	PURIMIX
		Cat. No.
PRODUCT	Main Set	CS-PMX-Set