

PURISPIN 6

Personal mini micro centrifuge



STANDARD SPECIFICATION

Description		PURISPIN 6
Max. RPM		6,200RPM (speed adjustments unavailable)
Max. RCF		1,800 xg
Max. Capacity		8 x 0.2/0.5/1.5/2.0ml tubes, 4 x 0.2ml PCR tube strips
Rotor	Angel Type	8 Holes Circular Fixed-Angel rotor for 0.5/1.5/2.0 Tubes
	Strip Type	32 Hole Strip Rotor for 0.2 PCR Strip or individual Tubes
Speed Control		Automatic/Constant speed
Driving Motor		DC motor
Permissible Temp. & Humidity		5°C~40°C, <80%
Rotating Direction		Clockwise
Adaper		0.2/0.5ml
Power Consumption		8W
Safety Device		Door Limit Switch
Metrial	Lid	PC (Polycarbonate)
	Body	ABS
	Rotor	ABS
Deimension (W x D x H) & Net Weight		155 x 183 x 125mm, 1.34Kg
Packing Size (W x D x H) & Gross Weight		190 x 250 x 215mm, 1.60Kg
Power Supply		1 Phase, AC 120V, 60Hz or AC 230V, 50/60 Hz
Cat. No.		CS-C-M6-SET

MAIN FEATURE



Compact Design



Excellent Durability



Simple Operation



Silent Operation

ACCESSORIES



Fixed Angle Rotor



Fixed Angle Rotor for PCR strip



Adaptors for 0.2ml and 0.4/0.5ml tubes

ORDER INFORMATION

ITEM	DESCRIPTIONS	Cat. No.
PURISPIN 6 Set	120/230V, 50/60Hz, 1Phase, EU Plug, 155x183x125mm, 1.34kg Including 0.2ml / 0.5ml tube adapter, PCR Strips	CS-C-M6-SET